

Schircks Laboratories
 Postfach
 CH-8494 Bauma
 Switzerland
 Telephone +41 (0) 55 / 212 23 24
 E-mail labschircks@gmail.com
 Website www.schircks.ch


Schircks Laboratories

DATA SHEET

(6R,S)-5,10-Methylene-5,6,7,8-tetrahydrofolic acid, calcium salt

Abbreviation used: MethyleneFH₄

Product no. 16.226



C₂₀H₂₁N₇O₆(Ca)

MW: 495.5

Description	Yellow to brown crystalline powder		
Biochemistry	When methionine intake is low, levels of S-adenosylmethionine decrease. This stimulates the enzyme methylene tetrahydrofolate reductase (MTHFR) to convert 5,10-methylene tetrahydrofolate to 5-methyl-tetrahydrofolate, which is required for methionine synthesis.		
Solubility	MethyleneFH ₄ calcium salt is very slightly soluble in water. The solubility of methyleneFH ₄ is approx. 0.05 g per 100 g of water (22°C). The pH of 1 mg/ml solution is 6.6. Ultrasonication may be used to improve dissolution.		
Analytical methods	HPLC conditions:	column:	Waters Spherisorb S5-ODS1
		eluant:	10 mM Na ₂ HPO ₄ pH 7 / Methanol (10%)
		flow rate:	1 ml/min
		wavelength:	254 nm
		solution:	1 mg/3 ml in 0.05 M NaOH
	TLC: MethyleneFH ₄ solutions are too unstable to perform TLC.		
Specifications	Purity: HPLC	approximately 93%	
Stability	MethyleneFH ₄ solutions are unstable in the presence of air. 1 mM solutions in water are degraded by about 5% within 60 minutes of preparation whereas 1 mM solutions in 0.01 M NaOH are only slightly degraded after 1 hour of preparation. The solutions are pH sensitive. Solutions of pH 8 or higher are more stable, than those of pH less than 8. So it is preferable to make solutions in evacuated, dilute NaOH, (0.05 M NaOH) and use them immediately. The powder is relatively stable in the presence of air. The powder is degraded by 4% stored open at room temperature for 1 month. MethyleneFH ₄ is hygroscopic. Despite its instability in solution, the powder can be stored for several years, dry, in tightly closed vials at -20°C or colder.		
Storage	Keep the powder in ampoules at -20°C or colder. After an ampoule has been opened, the remainder should be stored in a tightly closed vial, in a freezer. MethyleneFH ₄ in ampoules can be transported without the use of dry ice. Dry, in ampoules, it is stable for several weeks at room temperature.		
Uses	MethyleneFH ₄ is an important standard for analytical work. It is sold for laboratory use only.		
Safety information	MethyleneFH ₄ is known to be safe and there are no special precautions required in handling this product.		
References	5,10-Methylene-5,6,7,8-tetrahydrofolate. Conformation of the Tetrahydropyrazine and Imidazolidine Rings, Martin Poe et al., Biochemistry, <u>18</u> No. 24,(1979), 5527. Synthesis of (6R,11S)- and (6R,11R)-5,10-Methylene[11-1H,2H]tetrahydrofolate. Stereochemical Paths of Serine Hydroxymethyltransferase, 5,10-Methylenetetrahydrofolate Dehydrogenase, and Thymidylate Synthetase, Lawrence J. Sliker and Stephen J. Benkovic, Catalysis J. Am. Chem. Soc. <u>106</u> 1, (1984), 833.		

Further data sheets can be found on our website www.schircks.ch

The information in this publication is based on our current knowledge and experience. It does not relieve users or processors from carrying out their own precautions and tests.