

4-Amino-7,8-dihydro-L-biopterin

Product number 11.281

Step 1: Reduction of 4-amino-L-biopterin

1.5 g of 4-amino-L-biopterin and 6.0 g of Na₂S₂O₄ are placed in a 250 ml round bottom flask. 60 ml of water are added and a glass stopcock is mounted. The round bottom flask is evacuated and swirled for a short time. Attention, the suspension can foam.

The round bottom flask is warmed in a water bath at 60°C for 40 minutes. The solution becomes light yellow.

HPLC analysis shows that this solution still contains about 6% of 4-amino-L-biopterin. More Na₂S₂O₄ or longer reaction time did not reduce the amount of 4-amino-L-biopterin.

The round bottom flask is immersed in cold water for 20 minutes. If the purification cannot be started immediately, the round bottom flask has to be stored at 2°C overnight.

In contrast to dihydrobiopterin, 4-amino-7,8-dihydro-L-biopterin is readily soluble in water and must therefore be purified by column chromatography.

Step 2: Column chromatography (Whatman CC 31)

Preparation of the column:

Column diameter: 8 cm, height of the stationary phase: 20 cm

A slurry of cellulose CC31 from Whatman and NH₄AcOH solution 5% is poured through a funnel into the column. The top of the cellulose is protected by filter papers. The column is equilibrated with H₂O.

The solution with the crude 4-amino-7,8-dihydro-L-biopterin (60 ml) is loaded onto the column.

The column is run with H₂O.

The fractions containing Na₂S₂O₄ are recognized by the smell.

The fractions are examined by TLC.

TLC foil: cellulose 400, developing solvent: water

The spots can be visualized by using ultraviolet light.

Solutions of 4-amino-7,8-dihydro-L-biopterin are very sensitive to oxygen.

They must therefore be protected by argon or stored in a freezer. All work must be done as quickly as possible.

When a container for the next fraction is changed, a small sample (3 drops) should be placed in a small vial and this vial and the container with the fraction should be placed immediately in a freezer.

The TLC is performed with the small samples in the vials.

The best fractions are collected and the water is evaporated to about 30 ml.

Step 3: Column chromatography (Sephadex G-10)

Column diameter: 5 cm, height of the stationary phase: 20 cm

The column is washed with 0.2% NaOH and then with a lot of water.

The eluent (water or water containing 0.5% of AcOH) is degassed and argon is bubbled through the eluent from above.

Schircks Laboratories

The solution of step 2 (30 ml) is loaded onto the column.

In order to protect the 4-amino-7,8-dihydro-L-biopterin from oxygen, argon is bubbled through the fractions.

The fractions are collected as described in step 2 and tested by HPLC.

Column: Spherisorb S5 ODS1

Eluant: 10 mM Na₂HPO₄ pH 6.0 plus 15% methanol


Detection: 254 nm

The best fractions are collected, the water is evaporated and the residue dried in a vacuum desiccator over NaOH.

The columns have to be washed immediately, so that the remaining 4-amino-7,8-dihydro-L-biopterin and the impurities do not oxidize in the columns.

Purity: 97.5% (HPLC)

Description: light yellow



4-Amino-7,8-dihydro-L-biopterin

C₉H₁₄N₆O₂
238.24
238.117823
C 45.37% H 5.92% N 35.27% O 13.43%

Product no. 11.281

HPLC	Sample	2 mg/ml
	Column	Spherisorb S5 ODS 1
	Eluent	10 mM PB pH 6.0 plus 15% MeOH
	Flow	1 ml/min
	Detection	254 nm