

## 6-Hydroxymethylpterin

Product number 11.420

CAS number 712-29-8

**Attention:** *6-Formylpterin is highly sensitive to light! 1)  
6-Hydroxymethylpterin is less sensitive to light. 2)  
6-Formylpterin is not stable in basic solutions (Cannizzaro  
reaction?) 3)*

To 7.9 g of 6-formylpterin are added 1.4 l of water 1). The suspension is mixed and NaOH 50% (about 8.8 g) is added until all is dissolved. Immediately afterwards 3) HCl 1:1 is added until a slight precipitation remains. Through a dropping funnel 0.5 N HCl are added to bring the pH value to 9.0.

Reduction of the 6-formylpterin to 6-hydroxymethylpterin with NaBH<sub>4</sub>:

*With an excess of NaBH<sub>4</sub> 6-hydroxy-dihydropterin is formed. This happens probably only when the pH is lowered to 6. The oxidation products of 6-hydroxy-dihydropterin are not easily removed. Perhaps the amount of NaBH<sub>4</sub> can be reduced. Nitrogen is bubbled through the suspension to remove the hydrogen. Wear safety glasses!*

2.8 g of NaBH<sub>4</sub> are slowly added within about 5 hours in 3 portions. From time-to-time HCl 1:1 is added to keep the pH at about 9.0. The suspension is mixed overnight and nitrogen is furthermore bubbled through the suspension. The pH is lowered by the addition of HCl 1:1 to 6.0 and set aside for three hours. The supernatant is removed by suction and the remaining suspension with the precipitated 6-hydroxymethylpterin is filtered through an 8 cm filtering funnel. The filter cake is rinsed with 100 ml of water and then dried in a vacuum desiccator over NaOH to give 3.5 g of 6-hydroxymethylpterin.

Recrystallisation:

To the raw 6-hydroxymethylpterin (pulverized and sieved) 1.3 l of water are added. In order to dissolve the 6-hydroxymethylpterin NaOH 50% is added until pH 12.0 is reached.

1.0 g of activated carbon is added to the red solution and the mixture is stirred for 20 minutes.

The mixture is filtered through a 5 cm filtering funnel.

To the filtrate is added with stirring first HCl 1:1 and then slowly while stirring vigorously 0.5 N HCl until a pH of about 5.5 is reached.

In order not to get a too fine precipitate the solution is refluxed for 10 minutes and then is allowed to cool overnight.

The precipitated 6-hydroxymethylpterin is filtered through a 6 cm filtering funnel. The filter cake is slowly rinsed with 200 ml of water and dried in a vacuum desiccator over NaOH to give 2.5 g of 6-hydroxymethylpterin 2).

## Schircks Laboratories

Purity: 98.2% (HPLC)

Description: light yellow powder

**Data Sheet: There is a data sheet available for this compound.**

Data sheets can be found in the price list by clicking on the product number of your choice.